

Disposable Latex Exam Gloves

Features

- Beaded cuff for easy donning and prevention of roll back
- Excellent and tactile sense
- Textured for improved grip
- Low level of extractable proteins and chemical residue
- Food safe, suitable for handling all foodstuffs
- AQL 1.5



Quality Standards

- Complies with EN 455



Applications

 Medical Purpose / Examination

 Industrial purpose / PPE

 Laboratory

 Healthcare and nursing

 General housekeeping

 IT Industry



SKU	Size	Color	Package	Box Size	Carton Dimension
LGBEM1002*	XS-XL	White	100pcs/box,10boxes/ctn	230*125*60mm	315*258*245mm
LGBEM1006*	XS-XL	Light Yellow	100pcs/box,10boxes/ctn	230*125*60mm	315*258*245mm

※ The last "*" means Size:3-XS 4-S 5-M 6-L 7-XL.

